

## GEMÜ D41

### Pneumatically operated diaphragm valve



#### Features

- Suitable for use in hygienic and aseptic applications (CIP/SIP capable)
- Fast, safe and simple diaphragm and actuator mounting via a central gearbox due to EasyLock technology
- Low maintenance thanks to the encapsulated diaphragm, which requires no retightening
- Very high Kv values due to flow-optimised valve body
- Identification of the mounting angle (hash mark)
- Optical position indicator and transparent cap as standard
- Simple modular expansion possible with future-orientated automation components

#### Description

The pneumatically operated GEMÜ D41 diaphragm valve is designed for use in sterile applications.

The diaphragm hermetically separates the actuator from the working medium.

The GEMÜ D41 diaphragm valve with EasyLock technology is mounted entirely without loose components using a central gearbox. All actuator parts (except the seals and design elements) are made from stainless steel. The normally closed, normally open and double acting control functions are available. The valve has an optical position indicator with a transparent cap as standard.

#### Technical specifications

- **Media temperature :** -10 to 100 °C
- **Sterilization temperature:** max. 150 °C
- **Ambient temperature:** -10 to 80 °C
- **Operating pressure :** 0 to 10 bar
- **Nominal sizes:** DN 10 to 25
- **Body configurations:** Straight through body | T body
- **Connection types:** Clamp | Spigot
- **Connection standards :** ASME | DIN | EN | ISO
- **Body materials:** (316L), forging material | 1.4435 (316L), block material
- **Diaphragm materials:** PTFE/EPDM
- **Conformities:** ATEX | FDA | Reg. (EU) No. 10/2011 | Regulation (EC) No. 1935/2004 | USP

Technical data depends on the respective configuration




further information  
webcode: GW-D41



## Product description

### Construction



| Position | Name   | Materials                  |
|----------|--|----------------------------|
| 1        | Transparent cap  | PC                         |
| 2        | CONEXO actuator RFID chip<br>(see Conexo information)  |                            |
| 3        | Actuator   | 1.4301 / 1.4308            |
| 4        | Distance piece with leak detection hole                | 1.4308 / 1.4408            |
| 5        | Diaphragm  | PTFE/EPDM two-piece        |
| 6        | CONEXO diaphragm RFID chip<br>(see Conexo information) |                            |
| 7        | Valve body   | 1.4435, forged body (316L) |
| 8        | CONEXO body RFID chip<br>(see Conexo information)      |                            |
| 18       | Special stud bolts                                     |                            |
| 19       | Hexagon SW6  |                            |

## GEMÜ CONEXO

The interaction of valve components that are equipped with RFID chips and an associated IT infrastructure actively increase process reliability.



The diaphragm is read out using an RFID reader, the CONEXO pen. The CONEXO app, which can be installed on mobile devices, not only facilitates and improves the “installation qualification” process, but also makes the maintenance process much more transparent and easier to document. The app actively guides the maintenance technician through the maintenance schedule and directly provides them with all the information assigned to the valve, such as test reports, testing documentation and maintenance histories. The CONEXO portal acts as a central element, helping to collect, manage and process all data.

**You can find further information on GEMÜ CONEXO at:**

[www.gemu-group.com/conexo](http://www.gemu-group.com/conexo)

## Range overviews

### Range overview of surface finishes

Internal surface finishes for block material body <sup>1)</sup>

| Process contact surfaces | Mechanically machined <sup>2)</sup> |      | Electropolished            |      |
|--------------------------|-------------------------------------|------|----------------------------|------|
|                          | Hygiene class<br>DIN 11866          | Code | Hygiene class<br>DIN 11866 | Code |
| Ra ≤ 0.40 µm             | H4                                  | 1536 | HE4                        | 1537 |

| Process contact surfaces in accordance with ASME BPE <sup>3)</sup> | Mechanically machined <sup>2)</sup> |      | Electropolished              |      |
|--|-------------------------------------|------|------------------------------|------|
|  | ASME BPE surface designation        | Code | ASME BPE surface designation | Code |
| Ra max. = 0.51 µm<br>(20 µinch)                                    | SF1                                 | SF1  | -                            | -    |
| Ra max. = 0.38 µm<br>(15 µinch)                                    | -                                   | -    | SF4                          | SF4  |

Ra acc. to DIN EN ISO 4288 and ASME B46.1

- 1) Surface finishes of customized valve bodies may be limited in special cases.
- 2) Or any other finishing method that meets the Ra value (acc. to ASME BPE).
- 3) When using these surfaces, the bodies are marked according to the specifications of ASME BPE.  
The surfaces are only available for valve bodies which are made of materials (e.g. GEMÜ material codes 40, 41) and use connections (e.g. GEMÜ connection codes 59, 80, 88) according to ASME BPE.

### Range overview of valve bodies

#### Range overviews of spigots 17, 59, 60

| MG       | DN        | Connection type code |    |    |
|----------|-----------|----------------------|----|----|
|          |           | 17                   | 59 | 60 |
|          |           | Material code        |    |    |
|          |           | 40                   |    |    |
| <b>B</b> | <b>10</b> | X                    | -  | X  |
|          | <b>15</b> | X                    | X  | X  |
|          | <b>20</b> | X                    | X  | -  |
|          | <b>25</b> | -                    | X  | -  |

MG = diaphragm size, X = standard

## Order data

The order data provide an overview of standard configurations.

Please check the availability before ordering. Other configurations available on request.

## Order codes

| 1 Type   | Code |
|--|------|
| Diaphragm valve, pneumatically operated, stainless steel piston actuator, EasyLock | D41  |

| 2 DN, connection 1 | Code |
|--------------------|------|
| DN 10              | 10   |
| DN 15              | 15   |
| DN 20              | 20   |
| DN 25              | 25   |

| 3 Body configuration | Code |
|----------------------|------|
| 2/2-way body         | D    |

| 4 Valve body connection type, connection 1   | Code |
|--|------|
| Spigot EN 10357 series A/DIN 11866 series A  | 17   |
| Spigot ASME BPE/DIN EN 10357 series C (from 2022 edition)/DIN 11866 series C   | 59   |
| Spigot ISO 1127/DIN EN 10357 series C (2014 edition)/DIN 11866 series B  | 60   |
| Clamp ASME BPE, for pipe ASME BPE, face-to-face dimension FTF ASME BPE, length only for body configuration D                           | 80   |
| Clamp DIN 32676, for pipe DIN 11866 series B, face-to-face dimension FTF EN 558 series 7, length only for body configuration D         | 82   |
| Clamp ASME BPE, for pipe ASME BPE, face-to-face dimension FTF EN 558 series 7, length only for body configuration D                    | 88   |
| Clamp DIN 32676, for pipe DIN 11866 series A, face-to-face dimension FTF acc. to EN 558 series 7, length only for body configuration D | 8A   |
| Clamp DIN 32676, for pipe DIN 11866 series C, face-to-face dimension FTF ASME BPE, length only for body configuration D                | 8P   |
| Clamp DIN 32676, for pipe DIN 11866 series C, face-to-face dimension FTF EN 558 series 7, length only for body configuration D         | 8T   |

| 5 Valve body material       | Code |
|-----------------------------|------|
| 1.4435 (F316L), forged body | 40   |

| 6 Diaphragm material | Code |
|----------------------|------|
| PTFE/EPDM, two-piece | T1   |

| 7 Control function   | Code |
|----------------------|------|
| Normally closed (NC) | 1    |
| Normally open (NO)   | 2    |

| 7 Control function   | Code |
|--|------|
| Double acting (DA)   | 3    |
| Normally closed (NC), interface/hexagon 90° offset to the pipework direction | Q    |
| Double acting (DA), interface/hexagon 90° offset to the pipework direction   | T    |
| Normally open (NO), interface/hexagon 90° offset to the pipework direction   | U    |

| 8 Actuator spring set | Code |
|-----------------------|------|
| Standard spring set   | 1    |

| 9 Diaphragm size | Code |
|------------------|------|
| Diaphragm size B | B    |

| 10 Actuator size | Code |
|------------------|------|
| Actuator size 2  | 2    |

| 11 Valve body surface   | Code |
|---|------|
| Ra ≤ 0.4 µm for media wetted surfaces, in accordance with DIN 11866 H4, mechanically polished internal mechanically machined internal | 1536 |
| Ra ≤ 0.4 µm for media wetted surfaces, in accordance with DIN 11866 HE4, electropolished internal/external                            | 1537 |
| Ra max. 0.51 µm (20 µin.) for media wetted surfaces, in accordance with ASME BPE SF1, mechanically machined internal                  | SF1  |
| Ra max. 0.38 µm (15 µin.) for media wetted surfaces, in accordance with ASME BPE SF4, electropolished internal/external               | SF4  |

| 12 CONEXO   | Code |
|---|------|
| Integrated RFID chip for electronic identification and traceability | C    |
| Without   |      |

| 13 Special function | Code |
|---------------------|------|
| ATEX marking        | X    |

**Order example**

| Ordering option                            | Code | Description  |
|--|------|--|
| 1 Type                                     | D41  | Diaphragm valve, pneumatically operated, stainless steel piston actuator, EasyLock                         |
| 2 DN, connection 1                         | 15   | DN 15  |
| 3 Body configuration                       | D    | 2/2-way body   |
| 4 Valve body connection type, connection 1 | 17   | Spigot EN 10357 series A/DIN 11866 series A  |
| 5 Valve body material                      | 40   | 1.4435 (F316L), forged body  |
| 6 Diaphragm material                       | T1   | PTFE/EPDM, two-piece   |
| 7 Control function                         | 1    | Normally closed (NC)   |
| 8 Actuator spring set                      | 1    | Standard spring set  |
| 9 Diaphragm size                           | B    | Diaphragm size B   |
| 10 Actuator size                           | 2    | Actuator size 2  |
| 11 Valve body surface                      | 1537 | Ra ≤ 0.4 µm for media wetted surfaces, in accordance with DIN 11866 HE4, electropolished internal/external |
| 12 CONEXO                                  |      | Without  |
| 13 Special function                        | X    | ATEX marking   |

## Technical data

### Medium

**Working medium:** Corrosive, inert, gaseous and liquid media which have no negative impact on the physical and chemical properties of the body and diaphragm material.

**Control medium:** Inert gases


### Temperature

| Media temperature: | Diaphragm material  | Standard      |
|--------------------|---------------------|---------------|
|                    | PTFE/EPDM (code T1) | -10 to 100 °C |

**Ambient temperature:** -10 – 80 °C

**Control medium temperature:** 0 – 60 °C

**Sterilization temperature:** PTFE/EPDM (code T1): Max. 150 °C, maximum 180 min per cycle  
 The sterilization temperature is valid for steam (saturated steam) or superheated water (max. 6 bar). PTFE diaphragms can also be used as moisture barriers; however, this will reduce their service life. It is not possible to retighten the connecting screws between the actuator and the valve body. This also applies to PTFE diaphragms exposed to high temperature fluctuations. The maintenance cycles must be adapted accordingly. The following valve arrangement has proven itself for interfaces between steam and process pipes: A globe valve for shutting off steam pipes and a diaphragm valve as an interface to the process pipes.



**Storage temperature:** -10 – 40 °C

### Pressure

**Operating pressure:** 0 – 10 bar  
 Information on operating pressures applied on both sides and for high purity media on request.  
 The operating pressures apply at room temperature. In case of deviating temperatures, observe the pressure / temperature correlation.

#### Control function 1

| MG       | DN             | Actuator size | Diaphragm material |
|----------|----------------|---------------|--------------------|
|          |                |               | PTFE               |
| <b>B</b> | <b>10 - 25</b> | <b>2</b>      | 0-10.0             |

MG = diaphragm size

All pressures are gauge pressures. Operating pressure values were determined with static operating pressure applied on one side of a closed valve. Sealing at the valve seat and atmospheric sealing is ensured for the given values.

**Vacuum:** Can be used up to a vacuum of 70 mbar (absolute)

**Pressure rating:** PN 16

**Leakage rate:** Leakage rate A to P11/P12 EN 12266-1

**Filling volume:**

| Actuator size | Diaphragm size | Control function     |                      | Valve stroke |
|---------------|----------------|----------------------|----------------------|--------------|
|               |                | NC                   | NO                   |              |
| <b>2</b>      | <b>B</b>       | 0.06 dm <sup>3</sup> | 0.05 dm <sup>3</sup> | 7.5 mm       |

NC = filling volume in open position

**Kv values:**

| DN        | Diaphragm size | Connection type (code) |      |      |
|-----------|----------------|------------------------|------|------|
|           |                | 17                     | 59   | 60   |
| <b>10</b> | <b>B</b>       | 4.4                    | -    | tbd. |
| <b>15</b> |                | 9.5                    | 3.9  | 9.5  |
| <b>20</b> |                | 9.5                    | 9.5  | -    |
| <b>25</b> |                | -                      | 10.0 | -    |

MG = diaphragm size, Kv values in m<sup>3</sup>/h

Kv values determined in accordance with EN 60534, inlet pressure 5 bar, Δp 1 bar, stainless steel valve body and two-piece PTFE/EPDM diaphragm. The Kv values for other product configurations (e.g. other diaphragm or body materials) may differ. In general, all diaphragms are subject to the influences of pressure, temperature and process parameters, as a result of which the Kv values deviate beyond the tolerance limit of the standard. The Kv value curve (Kv value dependent on valve stroke) can vary depending on the diaphragm material and term of use.



**Control pressure:** 4.0 Up to 8.0 bar

## Product conformities

**Machinery Directive:** 2006/42/EC

**Pressure Equipment Directive:** 2014/68/EU

**Explosion protection:** ATEX (2014/34/EU), order code Special version X

**ATEX marking:** Gas:  II 2 G Ex h IIC T6 ... T3 Gb X  
Dust :  II -/2 D Ex h -/IIIC T150 °C -/Db X

**Food:** FDA  
Regulation (EC) No. 1935/2004 (only for material code C3, 40, 42, 41, 43)  
Regulation (EC) No. 10/2011  
USP Class VI

## Mechanical data

**Weight:**
**Body**

| MG       | DN        | Spigot               |      |      |
|----------|-----------|----------------------|------|------|
|          |           | Connection type code |      |      |
|          |           | 17                   | 59   | 60   |
| <b>B</b> | <b>10</b> | 0.28                 | -    | 0.28 |
|          | <b>15</b> | 0.28                 | 0.28 | 0.28 |
|          | <b>20</b> | 0.28                 | 0.28 | -    |
|          | <b>25</b> | -                    | 0.28 | -    |

Weights in kg  
MG = diaphragm size



**Weight:**
**Actuator (control function 1, NC)**


| MG       | Actuator size | Control function 1<br>(NC) | Control function 2<br>(NO) |
|----------|---------------|----------------------------|----------------------------|
| <b>B</b> | 2             | 1.78                       | 1.71                       |

Weights in kg  
MG = diaphragm size

**Angle of rotation:**


| MG       | DN        | Connection type (code) |       |       |
|----------|-----------|------------------------|-------|-------|
|          |           | 17                     | 59    | 60    |
|          |           | $\alpha$               |       |       |
| <b>B</b> | <b>10</b> | 36.3°                  | -     | 23.4° |
|          | <b>15</b> | 18.0°                  | 38.1° | 12.3° |
|          | <b>20</b> | 7.9°                   | 18.3° | -     |
|          | <b>25</b> | -                      | 3.0°  | -     |

MG = diaphragm size



Dimensions

Actuator dimensions




| Actuator size | MG | A     | ØB   | A1   | G     | M     |
|---------------|----|-------|------|------|-------|-------|
| 2             | B  | 163.0 | 65.0 | 42.0 | G 1/8 | M16x1 |

Dimensions in mm  
MG = diaphragm size  
\* CT = A + H1 (see body dimensions)

## Body dimensions

### Spigot DIN/EN/ISO/ASME (code 17, 59, 60), forged material (code 40)

Connection type spigot DIN/EN/ISO/ASME (codes 17, 59, 60)<sup>1)</sup>, forged material (code 40, 42)



| MG | DN      | Pipe standard |      |     | EN 10357 series A/DIN 11866 series A |      |     | ASME BPE/DIN 11866 series C |       |      | ISO 1127/EN 10357 series C/DIN 11866 series B |      |     |
|----|---------|---------------|------|-----|--------------------------------------|------|-----|-----------------------------|-------|------|---|------|-----|
|    |         |               |      |     |                                      |      |     |                             |       |      |   |      |     |
|    |         |               |      |     | Connection code                      |      |     |                             |       |      |   |      |     |
|    |         | 17            |      |     | 59                                   |      |     | 60                          |       |      |   |      |     |
| L  | c (min) | G1            | H1   | Ød  | s                                    | H1   | Ød  | s                           | H1    | Ød   | s   |      |     |
| B  | 10      | 108.0         | 25.0 | 0.8 | 7.3                                  | 13.0 | 1.5 | 5.57                        | 9.53  | 0.89 | 9.4   | 17.2 | 1.6 |
|    | 15      | 108.0         | 25.0 | 0.8 | 10.3                                 | 19.0 | 1.5 | 7.15                        | 12.70 | 1.65 | 11.45   | 21.3 | 1.6 |
|    | 20      | 108.0         | 25.0 | 0.8 | 12.3                                 | 23.0 | 1.5 | 10.33                       | 19.05 | 1.65 | -   | -    | -   |
|    | 25      | 108.0         | 25.0 | 0.8 | -                                    | -    | -   | 13.50                       | 25.40 | 1.65 | -   | -    | -   |

Dimensions in mm

MG = diaphragm size

#### 1) Connection type


Code 17: Spigot EN 10357 series A/DIN 11866 series A, formerly DIN 11850 series 2

Code 59: Spigot ASME BPE/DIN EN 10357 series C (from 2022 issue)/DIN 11866 series C

Code 60: Spigot ISO 1127/DIN EN 10357 series C (2014 issue)/DIN 11866 series B

### Clamp DIN/EN/ISO/ASME (codes 80, 82, 88, 8A, 8P, 8T), forged material (code 40)

Connection type clamp DIN/ASME (codes 80, 88, 8P, 8T) <sup>1)</sup>, forged material (code 40, 42)



| MG | DN | NPS  | ød1             |        | ød3             |        | H1   | L               |        |
|----|----|------|-----------------|--------|-----------------|--------|------|-----------------|--------|
|    |    |      | Connection type |        | Connection type |        |      | Connection type |        |
|    |    |      | 80, 8P          | 88, 8T | 80, 8P          | 88, 8T |      | 80, 8P          | 88, 8T |
| B  | 10 | 3/8" | -               | -      | -               | -      | -    | -               | -      |
|    | 15 | 1/2" | 9.40            | 9.40   | 25.0            | 25.0   | 7.2  | 88.9            | 108.0  |
|    | 20 | 3/4" | 15.75           | 15.75  | 25.0            | 25.0   | 10.3 | 101.6           | 117.0  |
|    | 25 | 1"   | 22.1            | 22.1   | 50.5            | 50.5   | 13.5 | 114.3           | 127.0  |

Dimensions in mm

MG = diaphragm size

#### 1) Connection type

Code 80: Clamp ASME BPE, face-to-face dimension FTF ASME BPE, length only for body configuration D

Code 88: Clamp ASME BPE, for pipe ASME BPE, face-to-face dimension FTF EN 558 series 7, length only for body configuration D

Code 8P: Clamp DIN 32676 series C, face-to-face dimension FTF ASME BPE, length only for body configuration D

Code 8T: Clamp DIN 32676 series C, face-to-face dimension FTF EN 558 series 7, length only for body configuration D

### Connection type clamp DIN/ISO (code 82) <sup>1)</sup>, forged material (code 40, 42)

| Connection type clamp DIN/ISO code 82 / forged material (code 10, 12) |    |      |                 |                 |      |                 |
|---|----|------|-----------------|-----------------|------|-----------------|
| MG  | DN | NPS  | ød1             | ød3             | H1   | L               |
|   |    |      | Connection type | Connection type |      | Connection type |
|   |    |      | 82              | 82              |      | 82              |
| B   | 10 | 3/8" | 14.0            | 25.0            | 9.4  | 108.0           |
|   | 15 | 1/2" | 18.1            | 50.5            | 11.4 | 108.0           |
|   | 20 | 3/4" | -               | -               | -    | -               |
|   | 25 | 1"   | -               | -               | -    | -               |

Dimensions in mm

MG = diaphragm size

#### 1) Connection type

Code 82: Clamp DIN 32676 series B, face-to-face dimension FTF EN 558 series 7, length only for body configuration D

**Connection type clamp DIN/EN (code 8A)<sup>1)</sup>, forged material (code 40, 42)**

| MG       | DN | NPS  | ød1             | ød3             | H1   | L               |
|----------|----|------|-----------------|-----------------|------|-----------------|
|          |    |      | Connection type | Connection type |      | Connection type |
|          |    |      | 8A              | 8A              |      | 8A              |
| <b>B</b> | 10 | 3/8" | 10.0            | 34.0            | 7.3  | 108.0           |
|          | 15 | 1/2" | 16.0            | 34.0            | 10.3 | 108.0           |
|          | 20 | 3/4" | 20.0            | 34.0            | 12.3 | 117.0           |
|          | 25 | 1"   | -               | -               | -    | -               |

Dimensions in mm

MG = diaphragm size

**1) Connection type**

Code 8A: Clamp DIN 32676 series A, face-to-face dimension FTF acc. to EN 558 series 7, length only for body configuration D

## Accessories



### **GEMÜ 12A0**

#### **Intelligent electrical position indicator**

Independent of the actuator size and control function, the GEMÜ 12A0 electrical position indicator, as an automation module, is compatible with all pneumatically operated process valves of the new valve generation and with quarter turn valves. Contactless position detection determines the valve position precisely, reliably and without being subject to wear. The current valve position is displayed via high visibility LEDs, and fed back via electrical signals. In addition to this, there is an integrated mechanical position indicator. Modern communication interfaces, an integrated sensor system and the GEMÜ app operating option are all features that characterize the innovative electrical position indicator.



### **GEMÜ 44A0**

#### **Multi-functional valve actuation**

Independent of the actuator size, the GEMÜ 44A0 multi-functional valve actuation, as an automation module, is compatible with all pneumatically operated process valves with single acting linear actuator of the new valve generation. Depending on the order variant and the set device functions, the connected process valves can be controlled conventionally open/closed (combi switchbox) or the valve position can be precisely controlled (positioner). Contactless position detection determines the valve position precisely, reliably and without being subject to wear. The current valve position is displayed via high visibility LEDs, and fed back via electrical signals. In addition to this, there is an integrated mechanical position indicator. Modern communication interfaces, an integrated sensor system and the GEMÜ app operating option are all features that characterize this innovative product.



GEMÜ Gebr. Müller Apparatebau GmbH & Co. KG  
Fritz-Müller-Straße 6-8, 74653 Ingelfingen-Criesbach, Germany  
Phone +49 (0) 7940 1230 · [info@gemue.de](mailto:info@gemue.de)  
[www.gemu-group.com](http://www.gemu-group.com)