

Certificate

Working medium down to minus 30° Celsius


Use of GEMÜ diaphragm valves of the types 650 up to diaphragm size 100, 673 and 653 up to diaphragm size 50 for working medium at -30° to 130° Celsius

We hereby confirm that the GEMÜ valve types listed above in the material versions listed below are suitable for use at media temperatures of -30° to 130° Celsius.

The certificate applies to the following material versions:

GEMÜ valve body material code 40/code 42 (forged body 1.4435), code 41/code 43 (block material 1.4435) and GEMÜ diaphragm 54/5M (PTFE/EPDM)

The service life depends on the application conditions, meaning that the maintenance intervals must be adjusted accordingly.



M. Barghoorn
Head of Global Technics

Ingelfingen, 16/07/2024